

# Dental Probes



1470



1471-1



1472-2



1473-3



1474-4



1475-5



1476-6



1477-7



1478-8



1479-9



1480-10




1481-11



1482-2  
Woodson

# Burnishers and Pocket Probes Scalers

## Polishing Instruments



**NEW  
LASER SYSTEM**



- New marking process.
- Laser system marking.
- No more paint peeling off.
- Permanent marking.
- Antiglare, low depth marking.
- Maximum tip resistance against breaking.

10-8-8-7-3-2-1



**1497**  
Fig CP-C

11-8-6-3



**1498**  
Fig CP-11

12-8-6-3



**1499**  
Fig CP-12

11.5-8-6-5-5-3-5



**1500**  
Fig CP-12S  
O.M.S

# Periodontal Pocket Probes

**NEW  
LASER SYSTEM**



**1501**  
Williams



**1502**  
Merrit A



**1503**  
Merrit B



**1504**  
Goldman  
Fox



**1505**  
University of  
Michigan



**1506**  
26 G